

# Multi-center analysis of CIREN occupant lumbar bone mineral density and correlation with age and fracture incidence

Ashley A. Weaver, Mona Saffarzadeh, Kristen M. Beavers,  
R. Caresse Hightower, Anna N. Miller, Joel D. Stitzel



Virginia Tech – Wake Forest University Center for Injury Biomechanics  
Wake Forest University School of Medicine, Winston-Salem, NC

## Center for Injury Biomechanics




# Crash Injury Research and Engineering Network

*CIREN The Nation's Largest Learning Laboratory*




**Occupant  
Demographics  
(i.e. Age)**


**Osteopenia/  
Osteoporosis  
(BMD)**



**Crash Characteristics**




**Occupant Contacts**



**Injury Causation Investigation**

# Last Year's Goals

- Lumbar BMD of 281 WFU occupants (109 M, 172 F) was quantified using phantom-less CT method
  - Of those classified as osteopenic (<145 mg/cc) using this technique
    - 64% are undiagnosed in CIREN
    - 36% are correctly classified in CIREN



*Weaver, Traffic Injury Prevention, AAAM, 2015.*

# Current Goals of Study

- Quantify lumbar BMD of 873 CIREN occupants (372 M, 501F) from 8 centers using phantom-less CT method
  - 873 CT scans obtained from CIREN database
    - 8 centers: Wake Forest, UVA/Fairfax, UAB, MCW, U of Maryland/Baltimore, Seattle, Michigan, San Diego
  - Additional CT scans requested from current CIREN centers




*Charles "McC." Mathias, Jr., National Study Center  
for Trauma and Emergency Medical Systems*

# BMD Analysis

## Phantom-less vBMD Analysis Workflow


Validate Phantom-less  
vBMD Technique

Pilot data:  
DXA & QCT Measures  
(50 subjects; 17M, 33F)



Quantify vBMD in  
CIREN Occupants

Collect vBMD & tissue  
measures for calibration  
to mg/cc




Correlate vBMD with  
CIREN Data

Relate vBMD to CIREN  
occupant demographics,  
comorbidities, injuries



# Lumbar Spine BMD Analysis of Pilot & CIREN CT Scans



# ROIs for BMD Calibration from HU to mg/cc

## Phantom calibration measures

Port 2: 0 mg/cc CaHA


Port 3: 75 mg/cc CaHA

Port 4: 150 mg/cc CaHA

## Fat-muscle calibration measures

Muscle: right psoas

Fat: subcutaneous, anterior



INTable Calibration Phantom

[http://www.image-analysis.com/intable\\_phantom](http://www.image-analysis.com/intable_phantom)




0 75 150  
mg/cc CaHA

# Phantom Calibration Method

Phantom calibration technique applied to each pilot subject scan:

Linear regression fit: [0, 75, 150] mg/cc CaHA vs Ports 2-4 ROI measures




# Fat-Muscle Calibration Method

Fat-muscle calibration technique applied to each pilot subject scan:


Assumed ground truth values: Fat = -69, Muscle = 77 mg/cc CaHA

Linear regression fit: [-69 77] mg/cc CaHA vs [Fat Muscle] ROI measures





# Calibration Method Comparison

- Pilot dataset shows Fat-Muscle Calibration  $\approx$  Phantom Calibration
- Fat-Muscle Calibration applied to phantom-less CIREN scans




# CT-based Predictive mg/cc Threshold for Osteopenia

- DXA/CT, 50 subject pilot data



# 2015 CIREN Study Population


- 873 occupants (372 M, 501 F) across 8 centers
- Abdominal CT
- Ages 16+



## CIREN Comorbidities

Diagnosed osteopenia/osteoporosis  
92 subjects

10.5%



781 subjects


89.5%

No documented osteopenia/osteoporosis

# CIREN Truth Table of Bone Quality: Documented Comorbidities vs CT-Predicted vBMD

$\geq 145$  mg/cc - Normal 

$< 145$  mg/cc - Osteopenia 




Normal Bone Quality

Osteopenia or Osteoporosis

Normal Bone Quality

Osteopenia or Osteoporosis



**CIREN Comorbidity**


# CIREN Subjects

## Documented Comorbidity in CIREN Database?

- No osteopenia/osteoporosis comorbidity
- Osteopenia/osteoporosis comorbidity


- Female
- + Male
- Female
- ⋯ Male

Mean Lumbar vBMD (mg/cc CaHA)



Age (years)



# Significantly Higher # Rib/Sternum Fractures in CIREN Occupants with vBMD < 145 mg/cc



# Greater proportion of occupants with <145 mg/cc BMD sustained thoracolumbar, lumbar, & thoracic vertebral body fx

■ No Fx – Normal,  $\geq 145$  mg/cc  
■ Fx - Normal,  $\geq 145$  mg/cc

■ No Fx - Osteopenic,  $< 145$  mg/cc  
■ Fx - Osteopenic,  $< 145$  mg/cc




# Summary & Conclusions

## Phantom-less CT BMD Estimation Technique

Osteopenia indicated for lumbar BMD < 145 mg/cc

- Associated with increased # rib/sternum fxs
- Associated with a greater incidence of thoracolumbar, thoracic, and lumbar spine fxs
- Of those classified as osteopenic using this technique
  - 60% are undiagnosed in CIREN
  - 40% are correctly classified in CIREN




Technique useful for osteopenia/osteoporosis classification of CIREN occupants & other prospective/retrospective BMD studies

## Center for Injury Biomechanics



# Ongoing Research

- Increase sample size further by collecting non-injury CTs that are not currently in the database from CIREN centers
- Investigate height, weight, and BMI correlation with lumbar vBMD in larger sample
  - Poor correlation in current data; contrasts positive correlations (BMD vs weight/BMI) reported in larger studies (*Felson, et al., 1993; Reid, 2002*)
- ***Goal: Measure lumbar vBMD in 1000+ CIREN occupants for correlation with age, fx incidence, height, weight, & crash conditions***



# Acknowledgments



Wake Forest School of Medicine Translational Science Institute  
Wake Forest University Translational Science Center

*Thank you to our committed collaborators!*



*Frank Pintar, Dale Halloway*



*Chris Michetti, Jeff Crandall, Thomas Hartka*



*Russell Griffin, Shannon Carroll*



*Charles "McC." Mathias, Jr., National Study Center  
for Trauma and Emergency Medical Systems*

*Tom Scalea, Pat Dischinger*

# Contact Info

Ashley Weaver, [asweaver@wakehealth.edu](mailto:asweaver@wakehealth.edu)

Joel Stitzel, [jstitzel@wakehealth.edu](mailto:jstitzel@wakehealth.edu)

Center for Injury Biomechanics

